	750W	850W	950W
	60s	50s	40s
	20 min – 2 hrs		

HOW TO USE

1. Remove gel pack from back support pocket and use a microwave to heat  or a refrigerator to cool  the gel pack according to the table shown
Do not exceed the recommended heating guidelines – the product should not be too hot to touch – Do Not Boil – Do Not Freeze
2. Use 10 second intervals to increase heat if necessary
3. When desired temperature is reached place the gel pack securely into the back support pocket
4. Repeat as required
5. Maximum overall heating time of 1 min 20s

CAUTION

- Do not apply gel pack directly to the skin
- People with known or suspected circulatory problems should not use this product except under the direct supervision of a doctor
- Remove the product if it becomes uncomfortable, if in doubt seek advice

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WARNING

- Only use the product as directed
- For external use only
- Do not puncture the gel pack
- If solution is swallowed contact your doctor immediately
- If the fluid comes into contact with the skin wash with water and soap and rinse thoroughly
- If the fluid comes into contact with the eyes, rinse opened eye for several minutes with running water – if symptoms persist consult your doctor
- Keep out of the reach of children
- Do not mix with combustible materials
- Replace gel pack after 1yr of use

STORAGE

- Keep away from extreme and direct heat
- Store in a dry place

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